



**13**  
**Section**

Plaster Instruments






194-010-22

Saw Complete 220 Volt  
1500 g. Weight

**Technical Data:**

Electrical Connection 220 Volt 50 Hz, Input 90 Watt,  
Protective Insulated, Anti- signal interference,  
Motor- safety Switch.



194-022-50  
194-022-65

Saw blade only for Plaster  
and synthetic bandages  
Ø50 mm For Thin Plaster  
Ø65 mm For Thick Plaster



194-025-11 Spanner only




194-110-20 Engel



Stile  
194-115-20 20 cm  
194-115-23 23 cm  
194-115-26 26 cm  
194-115-37 37 mm



194-112-18 Esmarch  
18,0 cm



194-122-35  
Bruns  
serrated

194-122-20 Bruns



194-130-21 Mod. Seutin  
21,0 cm



194-140-62 Stryker



194-143-45 Hennig



194-155-14 Knowles 14,0 cm



194-164-18 Wolff - Böhler 18,0 cm  
194-166-24 24,0 cm



194-170-15 Wolff-Böhler Miniatur



194-195-23 Lorenz  
23,0 cm



194-200-19 Lister  
Exentric 19,0 cm




194-215-20 20,0 cm



194-223-23 Bergmann  
23,0 cm



194-230-14 Lister  
14,0 cm  
194-230-18 18,0 cm  
194-230-20 20,0 cm



194-240-19 Smith Mod USA  
18,0 cm  
194-240-20 20,0 cm  
194-240-23 23,0 cm